

uniRAD solves the problem of multiple diagnostic viewers and disparate site log-ins experienced by many Radiologists and remote reading groups.

Integration of images and reports and seamless consolidation of data using eRAD's unique uniRAD solution provide a single point of access to information stored in numerous locations by different PACS and other information systems.

uniRAD automatically determines which configuration to use when communicating with a specific device, so the user enjoys a single, consistent user interface; a common work list; and convenient web accessibility.

uniRAD is comprehensive and efficient

Radiologists can:

- View images and reports from differing HL7 systems, such as RIS, EMR, HIS, and PM
- Access diverse DICOM devices, including parallel and legacy PACS archives
- Incorporate dictation and speech recognition systems

uniRAD single-viewer solution is cost effective

- Saves money by eliminating the need for multiple viewers, worklists and reporting systems
- Simplifies office management and workflow by consolidating data automatically
- Increases productivity by reducing the number of hours spent manually

"uniRAD provides our group a single presentation with a comprehensive patient view to read and report each study."


George Schutz M.D., Radiologist, Oklahoma

uniRAD is flexible and user configurable

While uniRAD accommodates various external systems, it does not mandate their existence. With embedded scheduling, dictation, transcription, report distribution and other capabilities, uniRAD provides much of the function found in other information systems. Depending on the implementation, these tools can co-exist with other solutions, allowing the quick incorporation of non-integrated imaging sites and temporary alliances or loosely associated users.

uniRAD offers:

- Multi-site, multi-vendor PACS integration
- Consolidated worklist for productivity enhancements
- Web accessibility to existing PACS installations without migration or replacement
- Comprehensive clinical review through automatic prior-image and report retrieval
- Automated data management and routing
- Multi-discipline diagnostic viewer
- Co-existent report dictation, speech recognition and transcription solutions
- Single-source, enterprise-wide data analysis and reporting
- Programmable interface for custom integration



When combined with eRAD PAC's proven Work list, Hub Server, and Viewer technology, along with its archiving services, uniRAD's feature list includes

- Web-based PACS
- Ubiquitous Access
- Follow-me user profiles
- Efficient data streaming
- 128-bit encrypted communications
- Automated Failover
- Programmable web-based interface support
- Off-site hosting
- Disaster recovery

The result is a single-enterprise diagnostic information-management solution with the potential to leverage existing investments in assets, increase productivity, and lower costs for reading groups of any size.

To learn more about how uniRAD can bring seamless connectivity to your organization and bring order and control to your business visit us at www.erad.com or Call: **864.234.7430**.